

MED12L antibody



Catalog Number: 112583

Product name

MED12L antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MED12L Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 260-280kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_053002). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112583(MED12L antibody) at dilution of 1:500