

# MED20 antibody



Catalog Number: 112587

## Product name

MED20 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

MED20 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HeLa cells. Observed molecular weight by Western-blot: 28kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of MED20 recombinant protein (Accession Number: NM\_004275). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112587(MED20 antibody) at dilution of 1:500

## Clonality