

# MED31 antibody



Catalog Number: 112591

## Product name

MED31 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

MED31 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human brain tissue. Observed molecular weight by Western-blot: 16kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of MED31 recombinant protein (Accession Number: NM\_016060). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112591(MED31 antibody) at dilution of 1:600

## Clonality