

MED4 antibody



Catalog Number: 112592

Product name

MED4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MED4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 34-36kd

Preparation

This antibody was obtained by immunization of MED4 recombinant protein (Accession Number: NM_014166). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112592(MED4 antibody) at dilution of 1:600