

MED6 antibody



Catalog Number: 112593

Product name

MED6 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MED6 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, A431 cells. Observed molecular weight by Western-blot: 31-33kd

Preparation

This antibody was obtained by immunization of MED6 recombinant protein (Accession Number: NM_005466). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112593(MED6 antibody) at dilution of 1:500