

# MIA3 antibody



Catalog Number: 112606

## Product name

MIA3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MIA3 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 250 kDa

## Preparation

This antibody was obtained by immunization of MIA3 recombinant protein (Accession Number: BC047116). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112606(MIA3 Antibody) at dilution of 1:1000



IP Result of anti-MIA3 (IP:Catalog No:112606, 5ug; Detection:Catalog No:112606 1:1000) with HeLa cells lysate 3000ug.