

MICALL2 antibody



Catalog Number: 112611

Product name

Polyclonal

MICALL2 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

MICALL2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 109 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of MICALL2 recombinant protein (Accession Number: NM_182924). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112611(MICALL2 Antibody) at dilution of 1:1000

