

# MICALL2 antibody



Catalog Number: 112611

---

## Product name

MICALL2 antibody

## Specificity

Human; other species not tested.

## Antibody description

MICALL2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 109 kDa

## Preparation

This antibody was obtained by immunization of MICALL2 recombinant protein (Accession Number: NM\_182924). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112611(MICALL2 Antibody) at dilution of 1:1000