

MEF2C antibody



Catalog Number: 112635

Product name

MEF2C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEF2C Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HeLa cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 63kd

Preparation

This antibody was obtained by immunization of MEF2C recombinant protein (Accession Number: BC026341). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112635(MEF2C antibody) at dilution of 1:400