

MEF2C antibody



Catalog Number: 112635

Product name

MEF2C antibody

Polyclonal

Ig Type

Specificity

Human, Mouse, Rat; other species not tested.

Rabbit IgG

Applications

Antibody description

MEF2C Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HeLa cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 63kd

ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of MEF2C recombinant protein (Accession Number: BC026341). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112635(MEF2C antibody) at dilution of 1:400

Clonality