

MEF2D antibody



Catalog Number: 112638

Product name

MEF2D antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEF2D Rabbit Polyclonal antibody. Positive WB detected in C2C12 cell, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of MEF2D recombinant protein (Accession Number: XM_017001314). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



C2C12 cell were subjected to SDS PAGE followed by western blot with Catalog No:112638(MEF2D Antibody) at dilution of 1:600



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:112638(MEF2D Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)