

MEK3 antibody



Catalog Number: 112644

Product name

MEK3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEK3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 37-40 kDa

Preparation

This antibody was obtained by immunization of MEK3 recombinant protein (Accession Number: NM_145109). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112644(MAP2K3 antibody) at dilution of 1:1000



IP Result of anti-MAP2K3 (IP:Catalog No:112644, 4ug; Detection:Catalog No:112644 1:1000) with HeLa cells lysate 3200ug.