

MEPE antibody



Catalog Number: 112658

Product name

MEPE antibody

Specificity

human; other species not tested.

Antibody description

MEPE Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A2780 cells, HeLa cells. Observed molecular weight by Western-blot: 43kd,98kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_020203). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112658(MEPE antibody) at dilution of 1:600