

MKI67IP antibody



Catalog Number: 112664

Product name

MKI67IP antibody

Specificity

human; other species not tested.

Antibody description

MKI67IP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of MKI67IP recombinant protein (Accession Number: NM_032390). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112664(MKI67IP antibody) at dilution of 1:500