

MKRN2 antibody



Catalog Number: 112669

Product name

MKRN2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MKRN2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A375 cells, HeLa cells. Observed molecular weight by Western-blot: 47kd

Preparation

This antibody was obtained by immunization of MKRN2 recombinant protein (Accession Number: NM_014160). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112669(MKRN2 Antibody) at dilution of 1:1000