

# MND1 antibody

Catalog Number: 112714

## Product name

MND1 antibody

## Immunogen

[Human MND1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

MND1 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in human liver tissue, human kidney tissue, human lung tissue. Observed molecular weight by Western-blot: 30-35 kDa

## Preparation

This antibody was obtained by immunization of MND1 recombinant protein (Accession Number: NM\_032117). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112714(MND1 antibody) at dilution of 1:500



IP Result of anti-MND1 (IP:Catalog No:112714, 4ug; Detection:Catalog No:112714 1:600) with mouse liver tissue lysate 6000ug.