

MNDA antibody



Catalog Number: 112715

Product name

MNDA antibody

Specificity

human; other species not tested.

Antibody description

MNDA Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of MNDA recombinant protein (Accession Number: NM_002432). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112715(MNDA antibody) at dilution of 1:1000