

# MNS1 antibody



Catalog Number: 112717

## Product name

MNS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MNS1 Rabbit Polyclonal antibody. Positive WB detected in rat testis tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 55-60 kDa

## Preparation

This antibody was obtained by immunization of MNS1 recombinant protein (Accession Number: NM\_018365). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112717(MNS1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:112717(MNS1 Antibody) at dilution of 1:50 (under 40x lens)