

# MOCS2 antibody

Catalog Number: 112725

## Product name

MOCS2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MOCS2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, human liver tissue, mouse heart tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of M<sup>o</sup>CS2 recombinant protein (Accession Number: NM\_004531). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112725(MocS2 antibody) at dilution of 1:200



Immunofluorescent analysis of HeLa cells using Catalog No:112725(MocS2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)