

MORF4L2 antibody



Catalog Number: 112735

Product name

MORF4L2 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

MORF4L2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human liver tissue, mouse liver tissue. Observed molecular weight by Western-blot: 30-32kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of MORF4L2 recombinant protein (Accession Number: NM_001142432). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112735(MORF4L2 antibody) at dilution of 1:800

Clonality