

MPP1 antibody



Catalog Number: 112747

Product name

MPP1 antibody

Immunogen

[Human MPP1 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

MPP1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of MPP1 recombinant protein (Accession Number: NM_002436). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112747(MPP1 antibody) at dilution of 1:600