

# MPP6 antibody



Catalog Number: 112750

---

## Product name

MPP6 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

MPP6 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 55-61kd

## Preparation

This antibody was obtained by immunization of MPP6 recombinant protein (Accession Number: XM\_017012316). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112750(MPP6 antibody) at dilution of 1:300