

MPV17L antibody



Catalog Number: 112758

Product name

MPV17L antibody

Specificity

human, mouse; other species not tested.

Antibody description

MPV17L Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells. Observed molecular weight by Western-blot: 10 kDa

Preparation

This antibody was obtained by immunization of MPV17L recombinant protein (Accession Number: NM_173803). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:112758(MPV17L antibody) at dilution of 1:1000