

# MTFR1 antibody



Catalog Number: 112773

## Product name

MTFR1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MTFR1 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, A549 cells. Observed molecular weight by Western-blot: 37kd

## Preparation

This antibody was obtained by immunization of MTFR1 recombinant protein (Accession Number: NM\_014637). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112773(MTFR1 antibody) at dilution of 1:1000



IP Result of anti-MTFR1 (IP:Catalog No:112773, 4ug; Detection:Catalog No:112773 1:800) with mouse testis tissue lysate 4800ug.