

MTMR6 antibody

Catalog Number: 112788

Product name

MTMR6 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MTMR6 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse heart tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 64-67 kDa, 72-76 kDa

Preparation

This antibody was obtained by immunization of MTMR6 recombinant protein (Accession Number: NM_004685). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112788(MTMR6 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:112788(MTMR6 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:112788(MTMR6

MTMR6 antibody



Catalog Number: 112788

antibody) at dilution of 1:100 (under 40x lens)