

MTR antibody

Catalog Number: 112795

Product name

MTR antibody

Immunogen

[Human MTR Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

MTR Rabbit Polyclonal antibody. Positive IF detected in A375 cells. Positive IHC detected in human pancreas tissue, human heart tissue. Positive WB detected in A375 cells, mouse pancreas tissue. Observed molecular weight by Western-blot: 140 kDa

Preparation

This antibody was obtained by immunization of MTR recombinant protein (Accession Number: NM_000254). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:112795(MTR Antibody) at dilution of 1:600




Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:112795(MTR Antibody) at dilution of 1:50 (under 10x lens)

MTR antibody

Catalog Number: 112795



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:112795(MTR Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of A375 cells using Catalog No:112795(MTR Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)