

MPZL3 antibody



Catalog Number: 112798

Product name

MPZL3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MPZL3 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, mouse brain tissue. Observed molecular weight by Western-blot: 70 kDa, 28 kDa

Preparation

This antibody was obtained by immunization of MPZL3 recombinant protein (Accession Number: NM_198275). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112798(MPZL3 Antibody) at dilution of 1:1000



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112798(MPZL3 Antibody) at dilution of 1:1000