

MRAP2 antibody



Catalog Number: 112800

Product name

MRAP2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MRAP2 Rabbit Polyclonal antibody. Positive IHC detected in human brain.

Preparation

This antibody was obtained by immunization of MRAP2 recombinant protein (Accession Number: NM_001346542). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human brain using Catalog No:112800(MRAP2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:112800(MRAP2 antibody) at dilution of 1:50 (under 40x lens)