

# MRFAP1 antibody



Catalog Number: 112803

## Product name

MRFAP1 antibody

## Immunogen

[Human MRFAP1 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

MRFAP1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 15kd

## Preparation

This antibody was obtained by immunization of MRFAP1 recombinant protein (Accession Number: NM\_001272053). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112803(MRFAP1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112803(MRFAP1 antibody) at dilution of 1:50 (under 10x lens)