

MRPL13 antibody



Catalog Number: 112812

Product name

MRPL13 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity


Human, Mouse, Rat, Monkey; other species not tested.

IHC: 1:20-1:200

Antibody description

MRPL13 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive IP detected in mouse lung tissue. Positive WB detected in HeLa cells, COS-7 cells, HL-60 cells, Jurkat cells, mouse lung tissue, rat lung tissue. Observed molecular weight by Western-blot: 23kd

Validations




Preparation

This antibody was obtained by immunization of MRPL13 recombinant protein (Accession Number: NM_014078). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112812(MRPL13 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:112812(MRPL13 antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:



MRPL13 antibody

Catalog Number: 112812

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:112812(MRPL13 antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-MRPL13 (IP:Catalog No:112812, 4ug; Detection:Catalog No:112812 1:500) with mouse lung tissue lysate 4000ug.