

# MRPP3 antibody



Catalog Number: 112838

## Product name

MRPP3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MRPP3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, L02 cells, mouse liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 65-67 kDa

## Preparation

This antibody was obtained by immunization of MRPP3 recombinant protein (Accession Number: XM\_017021835). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112838(KIAA0391 antibody) at dilution of



Immunofluorescent analysis of HepG2 cells, using KIAA0391 antibody Catalog No:112838 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).