

MRPP3 antibody



Catalog Number: 112838

Product name

MRPP3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MRPP3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, L02 cells, mouse liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 65-67 kDa

Preparation

This antibody was obtained by immunization of MRPP3 recombinant protein (Accession Number: XM_017021835). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112838(KIAA0391 antibody) at dilution of



Immunofluorescent analysis of HepG2 cells, using KIAA0391 antibody Catalog No:112838 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).