

MRPS31 antibody



Catalog Number: 112856

Product name

MRPS31 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MRPS31 Rabbit Polyclonal antibody. Positive IF detected in MCF-7. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 38kd

Preparation

This antibody was obtained by immunization of MRPS31 recombinant protein (Accession Number: NM_005830). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112856(MRPS31 antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells, using MRPS31 antibody Catalog No:112856 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).