

MRRF antibody



Catalog Number: 112862

Product name

MRRF antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MRRF Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive IF detected in MCF-7 cells. Positive WB detected in human lung tissue. Observed molecular weight by Western-blot: 29kd

Preparation

This antibody was obtained by immunization of MRRF recombinant protein (Accession Number: NM_138777). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112862(MRRF antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:112862(MRRF antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of MCF-7 cells, using

MRRF antibody



Catalog Number: 112862

MRRF antibody Catalog No:112862 at 1:50
dilution and Rhodamine-labeled goat anti-rabbit

IgG (red). Blue pseudocolor = DAPI (fluorescent
DNA dye).