

MTUS1 antibody



Catalog Number: 112883

Product name

MTUS1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MTUS1 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive IP detected in mouse brain tissue. Positive WB detected in Jurkat cells, human brain tissue, K-562 cells, PC-3 cells. Observed molecular weight by Western-blot: 40 kDa, 80 kDa, 110-120 kDa, 270-280 kDa

Preparation

This antibody was obtained by immunization of MTUS1 recombinant protein (Accession Number: BC033842). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112883(MTUS1 antibody) at dilution of 1:500



IP Result of anti-MTUS1 (IP:Catalog No:112883, 4ug; Detection:Catalog No:112883 1:500) with mouse brain tissue lysate 6000ug.



Immunofluorescent analysis of (-20oc Ethanol) fixed PC-3 cells using Catalog No:112883(MTUS1

MTUS1 antibody



Catalog Number: 112883

Antibody) at dilution of 1:50 and Alexa Fluor 488- conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)