

MUM1L1 antibody



Catalog Number: 112897

Product name

MUM1L1 antibody

Specificity

Human; other species not tested.

Antibody description

MUM1L1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human brain tissue.

Preparation

This antibody was obtained by immunization of MUM1L1 recombinant protein (Accession Number: BC109292). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:112897(MUM1L1 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:112897(MUM1L1 Antibody) at dilution of 1:50 (under 40x lens)