

Product name

MURC antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MURC Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 48-53 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001018116). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112899(MURC Antibody) at dilution of 1:600