

# MYCL1 antibody



Catalog Number: 112917

## Product name

Polyclonal

MYCL1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

MYCL1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Dilutions

## Preparation

Recommended Dilution:

This antibody was obtained by immunization of MYCL1 recombinant protein (Accession Number: BC011864). Purification method: Antigen affinity purified.

WB: 1:200-1:2000

## Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Validations

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112917(MYCL1 Antibody) at dilution of 1:700