

# MYEOV antibody



Catalog Number: 112924

## Product name

MYEOV antibody

## Specificity

Human; other species not tested.

## Antibody description

MYEOV Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in SGC-7901 cells, Jurkat cells. Observed molecular weight by Western-blot: 55-60 kDa

## Preparation

This antibody was obtained by immunization of MYEOV recombinant protein (Accession Number: XM\_047426763). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:112924(MYEOV antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:112924(MYEOV antibody) at dilution of 1:50 (under 10x lens)