

MYH3 antibody



Catalog Number: 112933

Product name

MYH3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MYH3 Rabbit Polyclonal antibody. Positive WB detected in mouse embryo tissue. Positive IP detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 220-230 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002470). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



mouse embryo tissue were subjected to SDS PAGE followed by western blot with Catalog No:112933(MYH3 Antibody) at dilution of 1:1000



IP Result of anti-MYH3 (IP:Catalog No:112933, 4ug; Detection:Catalog No:112933 1:300) with mouse skeletal muscle tissue lysate 4000ug.