

MYH6 antibody

Catalog Number: 112935

Product name

MYH6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYH6 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, mouse heart tissue. Positive WB detected in rat heart tissue. Observed molecular weight by Western-blot: 200-220 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002471). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IHC: 1:20-1:200


Validations



rat heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112935(MYH6 Antibody) at dilution of 1:5000



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112935(MYH6 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112935(MYH6 Antibody) at dilution of 1:50 (under 40x lens)