

NARG1L antibody



Catalog Number: 112954

Product name

NARG1L antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NARG1L Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 97-100kd

Preparation

This antibody was obtained by immunization of NARG1L recombinant protein (Accession Number: BC028112). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112954(NARG1L antibody) at dilution of 1:1000