

Nav1.5 antibody

Catalog Number: 112964

Product name

Nav1.5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Nav1.5 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 227 kDa

Preparation

This antibody was obtained by immunization of Nav1.5 recombinant protein (Accession Number: NM_001160160). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112964(SCN5A Antibody) at dilution of 1:600