

MYO1B antibody

Catalog Number: 112983

Product name

MYO1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYO1B Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, mouse brain tissue.

Preparation

This antibody was obtained by immunization of MYO1B recombinant protein (Accession Number: NM_001161819). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human heart using Catalog No:112983(MYO1B antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:112983(MYO1B antibody) at dilution of 1:50 (under 40x lens)