

# MYO1D antibody

Catalog Number: 112984

## Product name

MYO1D antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

MYO1D Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue, mouse brain tissue. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of MYO1D recombinant protein (Accession Number: BC030602). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112984(MYO1D antibody) at dilution of 1:200



IP Result of anti-MYO1D (IP:Catalog No:112984, 4ug; Detection:Catalog No:112984 1:300) with mouse lung tissue lysate 4000ug.