

# MYOM1-Specific antibody



Catalog Number: 112993

## Product name

MYOM1-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MYOM1-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse brain tissue, SKOV-3 cells. Positive IHC detected in human skeletal muscle tissue, human heart tissue. Observed molecular weight by Western-blot: 190kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003803). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112993(MYOM1-Specific antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:112993(MYOM1-Specific Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog

# MYOM1-Specific antibody



Catalog Number: 112993

---

No:112993(MYOM1-Specific Antibody) at dilution of 1:50 (under 40x lens)