

Myozenin 2 antibody



Catalog Number: 112999

Product name

Myozenin 2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Myozenin 2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 30-35kd

Preparation

This antibody was obtained by immunization of Myozenin 2 recombinant protein (Accession Number: NM_016599). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112999(MYOZ2 antibody) at dilution of 1:400