

# NCK2 antibody

Catalog Number: 113037

---

## Product name

NCK2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NCK2 Rabbit Polyclonal antibody. Positive WB detected in U-937 cells. Observed molecular weight by Western-blot: 43 kDa

## Preparation

This antibody was obtained by immunization of NCK2 recombinant protein (Accession Number: XM\_017005105). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:113037(NCK2 Antibody) at dilution of 1:1000