

# NCOA6 antibody

Catalog Number: 113044

## Product name

NCOA6 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NCOA6 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 68 kDa; 219 kDa; 110 kDa

## Preparation

This antibody was obtained by immunization of NCOA6 recombinant protein (Accession Number: NM\_001318240). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:500-1:5000

## Validations



IP Result of anti-NCOA6 (IP:Catalog No:113044, 4ug; Detection:Catalog No:113044 1:1000) with HEK-293 cells lysate 2560ug.