

NCOR2 antibody

Catalog Number: 113047

Product name

NCOR2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NCOR2 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in MCF7 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 180kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006312). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113047(NCOR2 antibody) at dilution of 1:300



Immunofluorescent analysis of HeLa cells using Catalog No:113047(NCOR2 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG